## Immunoprecipitation (1)

### Summary

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### FG beads® information

Product name	Protein G beads
Product number	TAS8848N1173
Storage temperature	2-8°C
Storage buffer	10mM HEPES(pH7.9), 50mM KCl, 1mM EDTA, 10%glycerol
Size of beads	190nm ± 20nm
IgG binding capacity	>100ug mouse IgG /mg of beads

### Method 1 (Binding Antibody)


1. **Wash**  
Transfer each amount of beads to a tube.  
Wash beads with 200ul PBS(-) 2 times at 4°C.
2. **Bind Antibody**  
Add antibody (5ug) diluted in 200ul PBS(-) to beads.
3. **Reaction**  
Mix for 30min at room temperature.
4. **Wash**  
Wash antibody binding beads with Wash buffer A 2 times at 4°C.

### Method 2 (Immunoprecipitation)


1. **Add sample solution**  
Add 200ul HeLa cell extracts containing HSA to the antibody binding beads.
2. **Reaction**  
Resuspend beads and incubate with rotation for 120min at 4°C.
3. **Wash**  
Separate magnetically and remove supernatant.  
Wash beads with Wash buffer B 3 times at 4°C.  
Transfer 100ul resuspended solution to new tube to perform 2 elution pattern.  
Separate magnetically and remove supernatant.
4. **Elution**
  - Acid elution**  
Add 28ul Elution buffer and resuspend beads.  
Stand for 5min at 4°C, separate magnetically.  
Transfer the supernatant to new clean tube and add 2ul Neutralization buffer, and 10ul 4 × SDS sample buffer.
  - Boil elution**  
Add 40ul 1 × SDS sample buffer and resuspend beads.  
Boil for 5min and remove the beads.
5. **Analyze the samples by SDS-PAGE and silver staining**

## Result (2)

### 1. Experiment with the same amount of beads

### 2. Experiment to obtain the same recovery amount of antigens

Acid elution samples show similar band pattern to Boil elution samples.