Immunoprecipitation (2)

Summary

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

Method 1 (Binding Antibody)

1. **Wash**
Transfer 0.1mg of beads to a tube.
Wash beads with 200ul PBS(-) 2 times at 4°C.
2. **Bind Antibody**
Add antibody (5ug) diluted in 200ul PBS(-) to beads.
3. **Reaction**
Mix for 30min at room temperature.
4. **Wash**
Wash antibody binding beads with Wash buffer A 2 times at 4°C.

Method 2 (Immunoprecipitation)

1. **Add sample solution**
Add 200ul HeLa cell extracts .
2. **Reaction**
Resuspend beads and incubate with rotation for 120min at 4°C.
3. **Wash**
Separate magnetically and remove supernatant.
Wash beads with Wash buffer B 3 times at 4°C.
Separate magnetically and remove supernatant.
4. **Elution**
Add 40ul SDS sample buffer and resuspend beads.
Boil for 5min and remove the beads.
5. **Analyze the samples by SDS-PAGE and silver staining**

FG beads® information

Product name	Protein G beads
Product number	TAS8848N1173
Storage temperature	2-8°C
Storage buffer	10mM HEPES(pH7.9), 50mM KCl, 1mM EDTA, 10%glycerol
Size of beads	190nm ± 20 nm
IgG binding capacity	>100ug mouse IgG /mg of beads

Methods 3 (Western blotting)

1. Perform SDS-PAGE, and place the gel in Transfer buffer for about 10min.
2. Transfer the protein from the gel to a PVDF membrane.
3. Block the membrane with the Blocking buffer for 15min at room temperature.
4. Dilute the primary antibody with the Blocking buffer to 1/150.
5. Incubate for 60min at room temperature.
6. Wash the membrane with TBS-T buffer 3 times.
7. Dilute the secondary antibody with the TBS-T buffer to 1/2000.
8. incubate for 60min at room temperature.
9. Wash the membrane with TBS-T buffer 5 times.
10. Detect with a chemiluminescence substrate.